


Proloquo4Text for Apple Watch Tutorial



Getting started with Proloquo4Text for Apple Watch



Requirements

Note that you need **watchOS 2** or later on your Apple Watch and **iOS 9** or later on your iPhone to have access to Proloquo4Text for Apple Watch, and **watchOS 3** to use text-to-speech output on your watch.

Set up your Apple Watch

To get started with your Apple Watch, pair and sync it with your iPhone. For more information, please see: <https://support.apple.com/en-us/HT204505>

- 1 Turn OFF Wake Screen on Wrist Raise**
Open the Settings app on Apple Watch, go to *General*, then *Wake Screen*, and turn OFF *Wake Screen on Wrist Raise* so Apple Watch will not dim when turning your wrist to show the message in Proloquo4Text to your communication partner.
- 2 Apple Watch language**
Make sure to set the language that you would like to use on Apple Watch. To do so, open Apple Watch companion app on the paired
- 3 iPhone**, tap *General*, tap *Language & Region*, and under the *System Language* section select your desired language

Configure Proloquo4Text for Apple Watch

In the Apple Watch companion app on your phone, **2** scroll down the list of installed apps to Proloquo4Text. Make sure **3** *Show App on Apple Watch* is turned ON. You can set whether text is displayed upside down or right side up **6** and whether the message will be displayed word by word or all at once **4**. When text is displayed word by word, you set the speed of the words in *Word's animation speed* **5**. Set whether the watch will speak aloud **7**, and which voice the watch will use. You can pick which of the Apple-supplied voices are used for each language. The Acapela voices used on your iPhone are not available on Apple Watch.


Displaying and speaking your message

When a sentence is selected it will be displayed full screen and upside-down so that it is easy to read for the communication partner when you flick your wrist in their direction. The message will be spoken automatically. To speak the message again, tap anywhere on the watch screen. To go back to form a new message, tap the back button on the upper left corner of the watch. **8**

Proloquo4Text for Apple Watch Tutorial



Customizing Proloquo4Text vocabulary for Apple Watch



Edit your Watch vocabulary

In the Proloquo4Text app, tap the name of name of the quick block ① until it says *Phrases*. Scroll down and tap the row labeled *Apple Watch*. ② Tap the pencil ③ to go into Edit mode.

④ Edit existing items

Tap any existing item to change its label or text to Speak.

⑤ Add items and categories

Tap Add, and then choose between ⑥ adding a new category for a group of items, or just a new item in the current category.

⑦ Syncing with your watch

When you're done with editing, tap *Done*. Your changes will automatically be sync'd to your watch.